

PD1-positive tumor-infiltrating lymphocytes (TIL) for the next generation of adoptive T cell therapy

Michelle R. Simpson-Abelson, Kenneth D'Arigo, Angel Cedano-Hilton, Ian Frank, Viktoria Gontcharova, Krit Ritthipichai, and Cécile Chartier

INTRODUCTION

Adoptive T cell therapy (ACT) with autologous tumor-infiltrating T cells (TIL) has demonstrated high response rates in patients with metastatic melanoma.¹ TIL products recognize tissue-specific antigens, neoantigens, and non-cancer related antigens.^{2,3} Neoantigen-specific T cells are considered the main contributors to the anti-tumor activity of TIL.⁴ Strategies enriching TIL for the neoantigen-specific T cells are expected to yield more potent therapeutic products, especially in epithelial cancers which contain a high proportion of non-cancer specific T cells.⁵ Several studies have demonstrated that expression of PD1 on TIL identifies the neoantigen-specific T cells.^{6,7,8} Presented here is the development of a new process to produce neoantigen-specific-enriched TIL products for clinical application.

MATERIALS & METHODS

- PD1-positive (PD1⁺) cells were sorted via flow cytometry directly from fresh tumor digests and expanded in vitro
- Samples from six melanomas, three sarcomas, six breast cancers, and eight lung cancers were evaluated
- 3 populations were studied:
 - PD1⁺ sorted TIL
 - PD1⁻ sorted TIL
 - Bulk TIL (whole tumor unsorted digest)
- TIL were evaluated for yield (cell count), phenotype (flow cytometry), TCRVβ repertoire (RNA-sequencing), non-specific functionality (anti-CD3 and PMA), and tumor reactivity and killing (co-culture assays)

Figure 1. A streamlined protocol for expanding PD1+ TIL to clinically relevant levels




Figure 1. The tumor is excised from the patient and transported to Iovance research laboratories. Upon arrival, the tumor is digested, and the single-cell suspension stained for CD3 and PD1. PD1+ TIL are sorted by FACS using an FX500 instrument (Sony). The PD1+ cell fraction is placed into a flask with an anti-human CD3 antibody (OCT3; 30ng/ml) and irradiated allogeneic PBMCs (feeders) at 1:100 (TIL:feeder) ratio and rapidly expanded for 22 days (REP).

- References**
- Rosenberg, S.A., et al., Durable complete responses in heavily pretreated patients with metastatic melanoma using T-cell transfer immunotherapy. *Clin Cancer Res*, 2011, 17(13): p. 4550-7.
 - Kvistborg, P., et al., TIL therapy broadens the tumor-reactive CD8(+) T cell compartment in melanoma patients. *Oncotarget*, 2012, 3(4): p. 409-418.
 - Simoni, Y., et al., Bystander CD8(+) T cells are abundant and phenotypically distinct in human tumor infiltrates. *Nature*, 2018, 557(7706): p. 575-579.
 - Schumacher, T.N. and R.D. Schreiber, Neoantigens in cancer immunotherapy. *Science*, 2015, 348(6230): p. 69-74.
 - Turcotte, S., et al., Phenotype and function of T cells infiltrating visceral metastases from gastrointestinal cancers and melanoma: implications for adoptive cell transfer therapy. *J Immunol*, 2013, 191(5): p. 2217-25.
 - Inozume, T., et al., Selection of CD8+PD-1+ lymphocytes in fresh human melanomas enriches for tumor-reactive T cells. *J Immunother*, 2010, 33(9): p. 956-64.
 - Gross, A., et al., PD-1 identifies the patient-specific CD8(+) tumor-reactive repertoire infiltrating human tumors. *J Clin Invest*, 2014, 124(5): p. 2246-59.
 - Thommen, D.S., et al., A transcriptionally and functionally distinct PD-1(+) CD8(+) T cell pool with predictive potential in non-small-cell lung cancer treated with PD-1 blockade. *Nat Med*, 2018.

DISCLOSURE & FUNDING STATEMENT

This study and poster are sponsored by Iovance Biotherapeutics, Inc. All authors are employees of Iovance Biotherapeutics, Inc.

RESULTS

Frequency of PD1+ TIL varies across tumor samples but in vitro expansion process reliably yields more than 1 billion TIL




Figure 2: Selected and bulk TIL were expanded from melanoma (n=6), lung cancer (n=7), breast cancer (n=6), and sarcoma (n=3). (A) Frequencies of PD1⁺ cells in fresh tumor digests are shown for each individual sample. Horizontal and vertical lines represent the mean values and standard errors, respectively. (B) PD1⁺ and PD1⁻ sorted cells, and bulk digests were expanded as described in Figure 1. Cells were counted at the completion of the REP and fold expansions (final cell count/seeding cell count) calculated that were used to extrapolate total cell counts. For Bulk TIL, seeding cell count was estimated using the percentage of T cells in the tumor digests. Mean values are plotted as bars and standard errors shown as vertical lines.

PD1+ TIL demonstrate a different phenotypic profile, compared to PD1- TIL







Figure 3. Digested tumors from melanoma (n=2), lung (n=2), and breast (n=2) were assessed phenotypically by flow cytometry, prior to sorting. (B-C) Live lymphocytes were gated on CD3⁺ cells and assessed for PD1⁺ and PD1⁻. The PD1⁺ and PD1⁻ populations were assessed for cell surface expression of (B) activation and (C) exhaustion markers. Mean values are plotted as bars and standard errors shown as vertical lines. Statistical significance was assessed by a paired student t-test: ****p<0.0001, *p<0.05.

In vitro expanded PD1+ TIL are phenotypically similar to bulk TIL



Expanded PD1+ TIL are oligoclonal and comprise a fraction of the clones present in bulk TIL



Expanded PD1+ TIL are functional as determined by IFNγ secretion and CD107a mobilization in response to non-specific stimulation




Figure 6: (A) PD1⁺-derived TIL, PD1⁻-derived TIL, and bulk TIL from melanoma (n=5), lung (n=6), and breast (n=6) were stimulated for 18 hours with plate-bound anti-CD3. Supernatants were assessed for IFN γ secretion by ELISA. Results are plotted for individual samples. (B) PD1⁺-derived TIL, PD1⁻-derived TIL, and bulk TIL from melanoma (n=5), lung (n=7), breast (n=6), and sarcoma (n=1) were assessed for CD107a cell surface expression in response to PMA stimulation for 4 hours on the CD4⁺ and CD8⁺ cells, by flow cytometry. Results are plotted for individual samples. Horizontal lines represent the mean percentages of each subset and vertical lines represent the standard errors.

Expanded PD1+ TIL demonstrate an enhancement in autologous melanoma cell killing and tumor reactivity relative to PD1- TIL




Figure 7: Tumor reactivity was assessed on PD1 selected TIL product from one melanoma sample. (A) Whole tumor digest was cleaned up using a dead cell removal kit (Miltenyi). 1e5 live cells were plated per well of a 96 well plate and permitted to adhere for 18 hours at 37°C in the xCELLigence instrument (ACEA Biosciences, Inc.). 1e5 PD1⁺- and PD1⁻-derived autologous TIL were added to their respective wells, resulting in a 1:1 (TIL:target) cell ratio, and incubated for 48 hours. Killing of the autologous target cells was recorded as increased impedance resulting from cell detachment. Cell killing (% cytotoxicity) (left most graph) was calculated using the formula $Cytotoxicity = [1 - (NCI_{t} / (AvgNCI_{t}))] \times 100$, where NCI_{t} is the Normalized Cell Index for the sample and NCI_{R} is the average of the Normalized Cell Index for the matching reference wells (digest alone). Right graph shows the normalized cell indices of the samples. (B) 1e5 cells from the whole tumor digest were cocultured with 1e5 TIL (or digest and TIL alone) for 18 hours. Supernatants were assessed for IFN γ release by ELISA (R&D systems). Bars represent the mean values of duplicate wells and vertical lines represent the standard errors.

SUMMARY

- Expanded PD1⁺ TIL demonstrate oligoclonality, compared to PD1⁻ derived TIL and bulk TIL, a sign of antigen-driven clonal expansion at the tumor site
- Preliminary data demonstrate autologous tumor cell killing by PD1⁺ but not PD1⁻ derived TIL
- Functionality of the expanded PD1⁺ TIL was confirmed by robust IFN γ and CD107a expression in response to non-specific stimulation
- Importantly, *in vitro* expansion of PD1⁺ TIL resulted in products phenotypically comparable with bulk TIL, indicating a strong therapeutic potential
- T cell markers regulated at the surface of expanded PD1⁺ TIL relative to pre-sort TIL included PD1 and CD25 and suggest a high activation level
- We intend to investigate this TIL product in clinic